


BPL XRAD₁₀₀

BPL X Rad 100 is a 100 mA, 100 pps Line Frequency Fixed X-ray solution that combines excellent Tube Head movement & intuitive User Interface and is ideal for Nursing Homes, Poly Clinics, Diagnostic Centres etc.

The Reliable Fixed X-Ray Imaging solution

- 1 User Friendly Design
- 2 Superior Image Quality
- 3 Maximum Value
- 4 Minimum Downtime




Powered By:


A O C

Advanced Operator Console

BPL X Rad 100 is powered by Advanced User friendly Operator Console with soft touch control of all radiographic parameters. The equipment uses State of the Art Microcontroller to control various functions and features.




- Alphanumeric LCD Display
- User Configurable Anatomical Programming
- 155 Anatomical Programs
- Pre-Set key for instant selection
- Self-Diagnostic Program
- X-Ray Counter with Password protection
- Memory to retain last selected factors
- 40 to 100 kVp in steps of 2 kVp only



Easy to use Counter Balanced Tube Head




Storage space for Accessories, Notepad etc.




Independent Voltmeter for monitoring input line Voltage



Independent Ammeter for indication of Tube Current



Anatomical Programming for precise selection of kV & mAs based on patient's anatomy



Self-Diagnostic Program to minimize downtime